

	PureVision <i>stoc permanent</i>	PureVision Therapeutice <i>stoc permanent</i>	SofLens Comfort <i>stoc permanent</i>	SofLens One Day <i>stoc permanent</i>
<i>Material</i>	Balafilcon A	Balafilcon A	Hilafilcon	Hilafilcon A
<i>Continutul de apa</i>	36%	36%	59.00%	70%
<i>Transmisibilitatea oxigenului Dk/t</i>	110	110	15.7	19
<i>Raza de curbura</i>	8.6 mm	8.6 mm	8.6 mm	8.6 mm
<i>Diametru</i>	14 mm	14 mm	14.2 mm	14 mm
<i>Ambalaj</i>	3 lentile	1 lentila	6 lentile	30 lentile
<i>Mod de purtare</i>	lunar (zi si noapte)	lunar (zi si noapte)	lunar	zilnic, de unica folosinta
<i>Dioptrii</i>	+ 6.00 D la - 6.00 D (pasul de 0.25D)	0.00 D	+ 6.00 D la - 6.00 D (pasul de 0.25D)	+ 6.50 D la - 6.00 D (pasul de 0.25D)
	- 6.50 D la - 12.00 D (pasul de 0.50 D)		- 6.50 D la - 9.00 D (pasul de 0.50 D)	- 6.50 D la - 9.00 D (pasul de 0.50 D)
	SofLens 66 Toric <i>comanda - 30 zile</i>	Natural Colors <i>stoc permanent (raza de 8.7)</i>	Optima FW <i>stoc permanent (raza de 8.7)</i>	Optima 38 <i>stoc permanent (raza de 8.7)</i>
<i>Material</i>	Alphafilcon A	Polymacon	Polymacon	Polymacon
<i>Continutul de apa</i>	66%	38.60%	38.60%	38.60%
<i>Transmisibilitatea oxigenului Dk/t</i>	16	14	24.3	14.2
<i>Raza de curbura</i>	8.5 mm	8.4, 8.7 mm	8.4, 8.7, 9.0 mm	8.4, 8.7, 9.0 mm
<i>Diametru</i>	14.5 mm	14 mm	14 mm	14 mm
<i>Ambalaj</i>	6 lentile	2 lentile plane / o lentila cu dioptrii	4 lentile	flacon cu 1 lentila
<i>Mod de purtare</i>	lunar	port ocazional (3 luni)	trimestrial (3 luni)	anual
<i>Dioptrii</i>	+ 6.00 D la - 6.00 D (pasul de 0.25D)	- 0.50 D la - 5.00 D (pasul de 0.25D)	+ 4.00 D la - 6.00 D (pasul de 0.25D)	+ 5.00 D la - 5.00 D (pasul de 0.25D)
	- 6.50 D la - 9.00 D (pasul de 0.50 D)	- 5.50 D la - 6.00 D (pasul de 0.50 D)	- 6.50 D la - 9.00 D (pasul de 0.50 D)	- 5.50 D la -12.00 D (pasul de 0.50 D)
<i>Cilindru</i>	- 0.75 D, -1.25 D, -1.75 D, - 2.25 D			
	de la 0.00 D la -9.00 D cilindru de - 2.75 D			
<i>Ax</i>	din 10° in 10°, de la 0° la 180°			