



# Bausch & Lomb

Perfecting Vision. Enhancing Life.™

## Arunca o privire: Intretinerea lentilelor de contact moi

### SOLUTII MULTIFUNCTIONALE

| Produs   | Descriere   | Dezinfectant | Surfactant  | Agent de deproteinizare | Mod de utilizare   | Mod de ambalare   |
|--|---|--------------|---|-------------------------|--|---|
| ReNu® MoistureLoc®<br>solutie multifunctionala | asigura confortul maxim zi de zi, curata, dezinfecteaza si indeparteaza depunerile de proteine, lipide si alte impuritati | Alexidina    | MoistureLoc™<br>(Poloxamer 407, Polyquaternium 10),<br>Poloxamina | Hydranate®              | se clateste fiecare fata a lentilei timp de minim 5 secunde, apoi lentilele se lasa in solutie minim 4 ore | flacon steril de plastic continand solutie la 60ml, 120ml, 240ml, 360ml |
| ReNu MultiPlus®<br>solutie multifunctionala    | curata, dezinfecteaza si indeparteaza depunerile de proteine, lipide si alte impuritati                                   | Dymed®       | Poloxamina  | Hydranate®              | se clateste fiecare fata a lentilei timp de minim 5 secunde, apoi lentilele se lasa in solutie minim 4 ore | flacon steril de plastic continand solutie la 60ml, 120ml, 240ml, 355ml |